

Linac Quality Assurance (TG-142)



THEUNIVERSITYofTENNESSEE



TG-142: Overview

The Report

- The Task Group 142 protocol was first published in 2009
- It is used for quality assurance of “C-Arm” linear accelerators
 - *Tomotherapy and Cyberknife have separate reports*
- It is used mechanical, dosimetry, and imaging testing
- The TG 142 report (2009) replaces the older TG40 QA report (1994)

Task Group 142 report: Quality assurance of medical accelerators^a

Eric E. Keall
Washington University, St. Louis, Missouri

Joseph Harkey
Rensselaer Polytechnic Institute Medical Center, Rensselaer, New York

John Bayouth
University of Iowa, Iowa City, Iowa

Felix Gershenson
Duke University, Durham, North Carolina

William Siron
University of Miami, Miami, Florida

Sean Dressler
Siteman Cancer Center, Atlanta, Georgia

Christopher Sengpiel
Mass. Gen. Hospital, Boston, Massachusetts

François Mignot
M. D. Anderson Cancer Center, Houston, Texas

Ulf Engström
University of California, San Francisco, San Francisco, California

Björn Arvidsson
KU Cancer Center, Houston, Texas

Chiharu Ito
Kumamoto University, Kumamoto, Japan

Consultants

Carlo Sandrin
Royal Marsden, London, United Kingdom

Todd Holmes
Varian Medical Systems, Palo Alto, California

Received 17 January 2009; revised 8 July 2009; accepted for publication 8 July 2009;

published online in *Physics in Medicine and Biology*, 17 August 2009.

The task group (TG) for quality assurance of medical accelerators was convened by the American Association of Medical Physics (AAMP) in 2002. The task group was formed by the Quality Assurance and Therapy Committees and the Quality Assurance and Oncology Improvement Subcommittee. The task group had two main charges. First is update, as needed, recommendations of the 1994 TG-40 report of the AAMP TG-40 on the quality assurance of electron beams, including recommendations for asymmetric jaws, epithelial collimation (EMC), and dynamic/beam wedges. The TG accomplished the update in TG-142. Second is to provide recommendations for quality assurance of medical accelerators that use intensity-modulated radiation therapy (IMRT) but also for imaging devices that are part of the linear accelerator. The report is intended to be a guide for the quality assurance of medical accelerators that deliver IMRT. The report was designed to account for the types of treatments delivered with the particular machine. The report is not intended to be a guide for quality assurance of other types of medical accelerators that use intensity-modulated radiation therapy (IMRT) require different test matrix structures. There are specific recommendations for IMRT quality assurance that are not included in this report. The report is intended to be a guide for action levels for the physicist to implement particular actions, whether they are inspection, scheduled action, or immediate and corrective action. The report is planned to be flexible for the physics department to use as a guide for quality assurance of medical accelerators. The report includes tables according to daily, monthly, and annual reviews, along with unique tables for wedge systems, EMC, and dynamic beam wedges. The report is intended to be a guide for the physicist to implement an QA program by the physicist in regards to building a QA team, establishing procedures, training of personnel, documentation, and end-to-end system checks. The submitted forms of this report have

1071 Med. Phys. 34 (9), September 2009 0094-2405(2009)09(9)1071-14\$35.00 © 2009 Am. Assoc. Phys. Med. 107




THEUNIVERSITYofTENNESSEE

TG-142: Overview

Background


- Several new technologies were developed between 1994 and 2009 not included in TG40



TG-142: Overview

Background



- Baseline dosimetric values entered into treatment planning system to characterize and/or model the treatment machine directly affect calculated plans
- Values can deviate from their baseline as a result of
 - Machine malfunction
 - Mechanical breakdown
 - Physical accidents
 - Component failure
 - Major component replacement
 - Gradual changes as a result of aging



TG-142: Overview

Background


- Baseline dosimetric values entered into treatment planning system to characterize and/or model the treatment machine directly affect calculated plans
- Values can deviate from their baseline as a result of
 - Machine malfunction
 - Mechanical breakdown
 - Physical accidents
 - Component failure
 - Major component replacement
 - Gradual changes as a result of aging



TG-142: Daily Tests

TABLE I. Daily.

Procedure	Machine-type tolerance		
	Non-IMRT	IMRT	SRS/SBRT
Dosimetry			
X-ray output constancy (all energies)			
Electron output constancy (weekly, except for machines with unique e-monitoring requiring daily)		3%	
Mechanical			
Laser localization	2 mm	1.5 mm	1 mm
Distance indicator (ODI) @ iso	2 mm	2 mm	2 mm
Collimator size indicator	2 mm	2 mm	1 mm
Safety			
Door interlock (beam off)		Functional	
Door closing safety		Functional	
Audiovisual monitor(s)		Functional	
Stereotactic interlocks (lockout)	NA	NA	Functional
Radiation area monitor (if used)		Functional	
Beam on indicator		Functional	



TG-142: Daily Tests

Dosimetry



X-ray output constancy (all energies)
 Electron output constancy (weekly,
 except for machines with unique
 e-monitoring requiring daily)

Mechanical

Laser localization
 Distance indicator (ODI) @ iso
 Collimator size indicator

Safety

Door interlock (beam off)
 Door closing safety
 Audiovisual monitor(s)
 Stereotactic interlocks (lockout)
 Radiation area monitor (if used)
 Beam on indicator



TG-142: Daily Tests

Dosimetry


X-ray output constancy (all energies)
 Electron output constancy (weekly,
 except for machines with unique
 e-monitoring requiring daily)

Mechanical

Laser localization
 Distance indicator (ODI) @ iso
 Collimator size indicator

Safety

Door interlock (beam off)
 Door closing safety
 Audiovisual monitor(s)
 Stereotactic interlocks (lockout)
 Radiation area monitor (if used)
 Beam on indicator



TG-142: Daily Tests

Dosimetry



X-ray output constancy (all energies)
Electron output constancy (weekly,
except for machines with unique
e-monitoring requiring daily)

Mechanical

Laser localization
Distance indicator (ODI) @ iso
Collimator size indicator

Safety

Door interlock (beam off)
Door closing safety
Audiovisual monitor(s)
Stereotactic interlocks (lockout)
Radiation area monitor (if used)
Beam on indicator



TG-142: Monthly Tests

TABLE II. Monthly.

Procedure	Machine-type tolerance		
	Non-IMRT	IMRT	SRS/SBRT
Dosimetry			
X-ray output constancy			
Electron output constancy		2%	
Backup monitor chamber constancy			
Typical dose rate ^a output constancy	NA	2% (@ IMRT dose rate)	2% (@ stereo dose rate, MU)
Photon beam profile constancy		1%	
Electron beam profile constancy		1%	
Electron beam energy constancy		2%/2 mm	





TG-142: Monthly Tests

Dosimetry

- TG-51 Calibration in Solid Water covers X-ray and Electron output consistency
- Backup monitor chamber constancy is read off the linac console during TG-51
- Dose rate consistency is measured during TG-51 by taking readings at the Maximum and half max dose rates

Dosimetry
X-ray output constancy
Electron output constancy
Backup monitor chamber constancy
Typical dose rate^a output constancy

Photon beam profile constancy
Electron beam profile constancy
Electron beam energy constancy




TG-142: Monthly Tests

Dosimetry

- Beam profile consistency can be measured with an array device, the EPID imager, a scanning water tank, or an single ionization chamber
- The measured profile is compared against a standard reference profile
- Electron beam energy consistency can be inferred from the profile shape

Dosimetry
X-ray output constancy
Electron output constancy
Backup monitor chamber constancy
Typical dose rate^a output constancy

Photon beam profile constancy
Electron beam profile constancy
Electron beam energy constancy




TG-142: Monthly Tests

Dosimetry

- Beam profile consistency can be measured with an array device, the EPID imager, a scanning water tank, or an single ionization chamber
- The measured profile is compared against a standard reference profile
- Electron beam energy consistency can be inferred from the profile shape

Dosimetry

X-ray output constancy
Electron output constancy
Backup monitor chamber constancy
Typical dose rate^a output constancy
Photon beam profile constancy
Electron beam profile constancy
Electron beam energy constancy



TG-142: Monthly Tests

TABLE II. Monthly.


Procedure	Non-IMRT	IMRT	SRS/SBRT
Mechanical			
Light/radiation field coincidence ^b		2 mm or 1% on a side	
Light/radiation field coincidence ^b (asymmetric)		1 mm or 1% on a side	
Distance check device for lasers compared with front pointer		1mm	
Gantry/collimator angle indicators (@ cardinal angles) (digital only)		1.0°	
Accessory trays (i.e., port film graticle tray)		2 mm	
Jaw position indicators (symmetric) ^c		2 mm	
Jaw position indicators (asymmetric) ^d		1 mm	
Cross-hair centering (walkout)		1 mm	
Treatment couch position indicators ^e	2 mm/1°	2 mm/1°	1 mm/0.5°
Wedge placement accuracy		2 mm	
Compensator placement accuracy ^f		1 mm	
Latching of wedges, blocking tray ^g		Functional	
Localizing lasers	±2 mm	±1 mm	< ± 1 mm



TG-142: Monthly Tests

Mechanical

- Light/radiation field coincidence^b
- Light/radiation field coincidence^c (asymmetric)
- Distance check device for lasers compared with front pointer
- Gantry/collimator angle indicators (@ cardinal angles) (digital only)
- Accessory trays (i.e., port film graticle tray)
- Jaw position indicators (symmetric)^d
- Jaw position indicators (asymmetric)^d
- Cross-hair centering (walkout)
- Treatment couch position indicators^e
- Wedge placement accuracy
- Compensator placement accuracy^f
- Latching of wedges, blocking tray^g
- Localizing lasers





THE UNIVERSITY OF TENNESSEE

TG-142: Monthly Tests

Mechanical

- Light/radiation field coincidence^b
- Light/radiation field coincidence^c (asymmetric)
- Distance check device for lasers compared with front pointer
- Gantry/collimator angle indicators (@ cardinal angles) (digital only)
- Accessory trays (i.e., port film graticle tray)
- Jaw position indicators (symmetric)**
- Jaw position indicators (asymmetric)^d
- Cross-hair centering (walkout)
- Treatment couch position indicators^e
- Wedge placement accuracy
- Compensator placement accuracy^f
- Latching of wedges, blocking tray^g
- Localizing lasers

Radiation Light Field 5x5



```

    Top-Right:2.3513 cm Delta:0.5133/mm.
    Bottom-Right:2.4352 cm Delta:-0.6480/mm.
    Right Side Distance:4.9865 cm Delta:-0.13472/mm.

    Left-Top:2.4854 cm Delta:-0.1456/mm.
    Right-Top:2.4904 cm Delta:-0.09635/mm.
    Top Distance:4.9758 cm Delta:-0.24191/mm.

    Left-Bottom:2.4794 cm Delta:0.20582/mm.
    Right-Bottom:2.4955 cm Delta:-0.044824/mm.
    Bottom Distance:4.9749 cm Delta:-0.25065/mm.

    Total Radiation-Light Field Misalignment:0.58435
  
```




THE UNIVERSITY OF TENNESSEE


TG-142: Monthly Tests

Mechanical

- Light/radiation field coincidence^a**
- Light/radiation field coincidence^b (asymmetric)
- Distance check device for lasers compared with front pointer
- Gantry/collimator angle indicators (@ cardinal angles) (digital only)
- Accessory trays (i.e., port film graticle tray)
- Jaw position indicators (symmetric)**
- Jaw position indicators (asymmetric)^d
- Cross-hair centering (walkout)
- Treatment couch position indicators^e
- Wedge placement accuracy
- Compensator placement accuracy^f
- Latching of wedges, blocking tray^g
- Localizing lasers



Radiation Light Field 10x10




THE UNIVERSITY OF TENNESSEE

TG-142: Monthly Tests

Mechanical

- Light/radiation field coincidence^a**
- Light/radiation field coincidence^b (asymmetric)
- Distance check device for lasers compared with front pointer
- Gantry/collimator angle indicators (@ cardinal angles) (digital only)
- Accessory trays (i.e., port film graticle tray)
- Jaw position indicators (symmetric)**
- Jaw position indicators (asymmetric)^d
- Cross-hair centering (walkout)
- Treatment couch position indicators^e
- Wedge placement accuracy
- Compensator placement accuracy^f
- Latching of wedges, blocking tray^g
- Localizing lasers



Radiation Light Field 15x15



THE UNIVERSITY OF TENNESSEE

TG-142: Monthly Tests

Mechanical

- Light/radiation field coincidence^b
- Light/radiation field coincidence^b (asymmetric)
- Distance check device for lasers compared with front pointer
- Gantry/collimator angle indicators (@ cardinal angles) (digital only)
- Accessory trays (i.e., port film graticle tray)
- Jaw position indicators (symmetric)^c
- Jaw position indicators (asymmetric)^d
- Cross-hair centering (walkout)
- Treatment couch position indicators^e
- Wedge placement accuracy
- Compensator placement accuracy^f
- Latching of wedges, blocking tray^g
- Localizing lasers




 THOMPSON CANCER
SURVIVAL CENTER

THE UNIVERSITY OF TENNESSEE | 

TG-142: Monthly Tests

Mechanical

- Light/radiation field coincidence^b
- Light/radiation field coincidence^b (asymmetric)
- Distance check device for lasers compared with front pointer
- Gantry/collimator angle indicators (@ cardinal angles) (digital only)
- Accessory trays (i.e., port film graticle tray)
- Jaw position indicators (symmetric)^c
- Jaw position indicators (asymmetric)^d
- Cross-hair centering (walkout)
- Treatment couch position indicators^e
- Wedge placement accuracy
- Compensator placement accuracy^f
- Latching of wedges, blocking tray^g
- Localizing lasers





 THOMPSON CANCER
SURVIVAL CENTER

THE UNIVERSITY OF TENNESSEE | 

TG-142: Monthly Tests

Mechanical

- Light/radiation field coincidence^b
- Light/radiation field coincidence^b (asymmetric)
- Distance check device for lasers compared with front pointer
- Gantry/collimator angle indicators (@ cardinal angles) (digital only)
- Accessory trays (i.e., port film graticle tray)**
- Jaw position indicators (symmetric)^c
- Jaw position indicators (asymmetric)^d
- Cross-hair centering (walkout)
- Treatment couch position indicators^e
- Wedge placement accuracy
- Compensator placement accuracy^f
- Latching of wedges, blocking tray^g
- Localizing lasers




THE UNIVERSITY OF TENNESSEE

TG-142: Monthly Tests

Mechanical

- Light/radiation field coincidence^b
- Light/radiation field coincidence^b (asymmetric)
- Distance check device for lasers compared with front pointer
- Gantry/collimator angle indicators (@ cardinal angles) (digital only)
- Accessory trays (i.e., port film graticle tray)
- Jaw position indicators (symmetric)^c
- Jaw position indicators (asymmetric)^d
- Cross-hair centering (walkout)**
- Treatment couch position indicators^e
- Wedge placement accuracy
- Compensator placement accuracy^f
- Latching of wedges, blocking tray^g
- Localizing lasers

Collimator Walk Out




THE UNIVERSITY OF TENNESSEE

TG-142: Monthly Tests

Mechanical

Light/radiation field coincidence^b
Light/radiation field coincidence^b (asymmetric)
Distance check device for lasers compared with
front pointer
Gantry/collimator angle indicators
(@ cardinal angles) (digital only)
Accessory trays (i.e., port film graticle tray)
Jaw position indicators (symmetric)^c
Jaw position indicators (asymmetric)^d
Cross-hair centering (walkout)
Treatment couch position indicators^e
Wedge placement accuracy
Compensator placement accuracy^f
Latching of wedges, blocking tray^g
Localizing lasers

Wedge Placement




TG-142: Monthly Tests

Mechanical

Light/radiation field coincidence^b
Light/radiation field coincidence^b (asymmetric)
Distance check device for lasers compared with
front pointer
Gantry/collimator angle indicators
(@ cardinal angles) (digital only)
Accessory trays (i.e., port film graticle tray)
Jaw position indicators (symmetric)^c
Jaw position indicators (asymmetric)^d
Cross-hair centering (walkout)
Treatment couch position indicators^e
Wedge placement accuracy
Compensator placement accuracy^f
Latching of wedges, blocking tray^g
Localizing lasers


Wedge Placement



TG-142: Monthly Tests

Mechanical

- Light/radiation field coincidence^b
- Light/radiation field coincidence^b (asymmetric)
- Distance check device for lasers compared with front pointer
- Gantry/collimator angle indicators (@ cardinal angles) (digital only)
- Accessory trays (i.e., port film graticle tray)
- Jaw position indicators (symmetric)^c
- Jaw position indicators (asymmetric)^d
- Cross-hair centering (walkout)
- Treatment couch position indicators^e
- Wedge placement accuracy
- Compensator placement accuracy^f
- Latching of wedges, blocking tray^g
- Localizing lasers



 THOMPSON CANCER
SURVIVAL CENTER

THE UNIVERSITY OF TENNESSEE 

TG-142: Monthly Tests

TABLE II. Monthly.

Procedure	Non-IMRT	IMRT	SRS/SBRT
Safety			
Laser guard-interlock test		Functional	
Respiratory gating			
Beam output constancy		2%	
Phase, amplitude beam control		Functional	
In-room respiratory monitoring system		Functional	
Gating interlock		Functional	

 THOMPSON CANCER
SURVIVAL CENTER

THE UNIVERSITY OF TENNESSEE 

TG-142: Monthly Tests

Safety

Laser guard-interlock test

Respiratory gating

Beam output constancy

Phase, amplitude beam control
In-room respiratory monitoring system
Gating interlock



THE UNIVERSITY OF TENNESSEE | OR

TG-142: Monthly Tests

Safety

Laser guard-interlock test

Respiratory gating

Beam output constancy

Phase, amplitude beam control
In-room respiratory monitoring system
Gating interlock



THE UNIVERSITY OF TENNESSEE | OR

TG-142: Monthly Tests

Safety

Laser guard-interlock test




Respiratory gating

Beam output constancy

Phase, amplitude beam control

In-room respiratory monitoring system

Gating interlock



THOMPSON CANCER
SURVIVAL CENTER

THE UNIVERSITY OF TENNESSEE |

TG-142: Annual Tests

TABLE III. Annual.

Procedure	Machine-type tolerance		
	Non-IMRT	IMRT	SRS/SBRT
Dosimetry			
X-ray flatness change from baseline		1%	
X-ray symmetry change from baseline		±1%	
Electron flatness change from baseline		1%	
Electron symmetry change from baseline		±1%	
SRS arc rotation mode (range: 0.5–10 MU/deg)	NA	NA	Monitor units set vs delivered: 1.0 MU or 2% (whichever is greater) Gantry arc set vs delivered: 1.0° or 2% (whichever is greater)




THOMPSON CANCER
SURVIVAL CENTER

THE UNIVERSITY OF TENNESSEE |

TG-142: Annual Tests

Dosimetry


X-ray flatness change from baseline
 X-ray symmetry change from baseline
 Electron flatness change from baseline
 Electron symmetry change from baseline
 SRS arc rotation mode
 (range: 0.5–10 MU/deg)



TG-142: Annual Tests

Dosimetry

X-ray flatness change from baseline
 X-ray symmetry change from baseline
 Electron flatness change from baseline
 Electron symmetry change from baseline
 SRS arc rotation mode
 (range: 0.5–10 MU/deg)



The beam flatness F is assessed by finding the maximum D_{\max} and minimum D_{\min} dose point values on the beam profile within the central 80% of the beam width and then using the relationship:


$$F = 100 \times \frac{D_{\max} - D_{\min}}{D_{\max} + D_{\min}}$$



TG-142: Annual Tests

Dosimetry

X-ray flatness change from baseline
 X-ray symmetry change from baseline
 Electron flatness change from baseline
 Electron symmetry change from baseline
 SRS arc rotation mode
 (range: 0.5–10 MU/deg)



The beam symmetry S is usually determined at z_{\max} , which represents the most sensitive depth for assessment of this beam uniformity parameter. A typical symmetry specification is that any two dose points on a beam profile, equidistant from the central axis point, are within 2% of each other. Alternately, areas under the z_{\max} beam profile on each side (left and right) of the central axis extending to the 50% dose level (normalized to 100% at the central axis point) are determined and S is then calculated from:

$$S = 100 \times \frac{\text{area}_{\text{left}} - \text{area}_{\text{right}}}{\text{area}_{\text{left}} + \text{area}_{\text{right}}}$$



TG-142: Annual Tests

TABLE III. Annual.

Procedure	Machine-type tolerance		
	Non-IMRT	IMRT	SRS/SBRT
X-ray/electron output calibration (TG-51)		±1% (absolute)	
Spot check of field size dependent output factors for x ray (two or more FSs)		2% for field size $<4 \times 4 \text{ cm}^2$, 1% $\geq 4 \times 4 \text{ cm}^2$	
Output factors for electron applicators (spot check of one applicator/energy)		±2% from baseline	
X-ray beam quality (PDD ₁₀ or TMR ₁₀ ²⁰)		±1% from baseline	
Electron beam quality (R_{50})		±1 mm	
Physical wedge transmission factor constancy		±2%	
X-ray monitor unit linearity (output constancy)	±2% $\geq 5 \text{ MU}$	±5% (2–4 MU), ±2% $\geq 5 \text{ MU}$	±5% (2–4 MU), ±2% $\geq 5 \text{ MU}$
Electron monitor unit linearity (output constancy)		±2% $\geq 5 \text{ MU}$	
X-ray output constancy vs dose rate		±2% from baseline	
X-ray output constancy vs gantry angle		±1% from baseline	
Electron output constancy vs gantry angle		±1% from baseline	
Electron and x-ray off-axis factor constancy vs gantry angle		±1% from baseline	
Arc mode (expected MU, degrees)		±1% from baseline	
TBI/TSET mode		Functional	
PDD or TMR and OAF constancy		1% (TBI) or 1 mm PDD shift (TSET) from baseline	
TBI/TSET output calibration		2% from baseline	
TBI/TSET accessories		2% from baseline	



TG-142: Annual Tests

X-ray/electron output calibration (TG-51)

Spot check of field size dependent
output factors for x ray
(two or more FSS)
Output factors for electron applicators
(spot check of one applicator/energy)
X-ray beam quality (PDD_{10} or TMR_{10}^{20})
Electron beam quality (R_{50})
Physical wedge transmission
factor constancy
X-ray monitor unit linearity
(output constancy)
Electron monitor unit linearity
(output constancy)
X-ray output constancy vs dose rate
X-ray output constancy vs gantry angle
Electron output constancy vs
gantry angle
Electron and x-ray off-axis factor
constancy vs gantry angle
Arc mode (expected MU, degrees)
TBI/TSET mode
 PDD or TMR and OAF constancy
TBI/TSET output calibration
TBI/TSET accessories



- Must be performed in liquid water using calibrated chambers as described in TG-51 report



TG-142: Annual Tests

X-ray/electron output calibration (TG-51)

Spot check of field size dependent
output factors for x ray
(two or more FSS)
Output factors for electron applicators
(spot check of one applicator/energy)
X-ray beam quality (PDD_{10} or TMR_{10}^{20})
Electron beam quality (R_{50})
Physical wedge transmission
factor constancy
X-ray monitor unit linearity
(output constancy)
Electron monitor unit linearity
(output constancy)
X-ray output constancy vs dose rate
X-ray output constancy vs gantry angle
Electron output constancy vs
gantry angle
Electron and x-ray off-axis factor
constancy vs gantry angle
Arc mode (expected MU, degrees)
TBI/TSET mode
 PDD or TMR and OAF constancy
TBI/TSET output calibration
TBI/TSET accessories



- Doesn't have to be measured in water, but you already have your tank setup
- Measure 5x5, 10x10, 15x15, and 20x20 for each photon mean energy (*typically 2-3 energies*)
- Measure 10x10 and 15x15 cones for each electron energy (*typically 5 energies*)
- Compare to baseline



TG-142: Annual Tests

X-ray/electron output calibration (TG-51)
Spot check of field size dependent output factors for x ray (two or more FSS)
Output factors for electron applicators (spot check of one applicator/energy)
X-ray beam quality (PDD₁₀ or TMR₁₀)
Electron beam quality (R_{50})
Physical wedge transmission factor constancy
X-ray monitor unit linearity (output constancy)
Electron monitor unit linearity (output constancy)
X-ray output constancy vs dose rate
X-ray output constancy vs gantry angle
Electron output constancy vs gantry angle
Electron and x-ray off-axis factor constancy vs gantry angle
Arc mode (expected MU, degrees)
TBI/TSET mode
PDD or TMR and OAF constancy
TBI/TSET output calibration
TBI/TSET accessories



- Doesn't have to be measured in water, but you already have your tank setup
- Measure 10x10 field size for open, 15, 30, 45, and 60 degree wedges for each photon beam energy
- Compare to baseline



TG-142: Annual Tests

X-ray/electron output calibration (TG-51)
Spot check of field size dependent output factors for x ray (two or more FSS)
Output factors for electron applicators (spot check of one applicator/energy)
X-ray beam quality (PDD₁₀ or TMR₁₀)
Electron beam quality (R_{50})
Physical wedge transmission factor constancy
X-ray monitor unit linearity (output constancy)
Electron monitor unit linearity (output constancy)
X-ray output constancy vs dose rate
X-ray output constancy vs gantry angle
Electron output constancy vs gantry angle
Electron and x-ray off-axis factor constancy vs gantry angle
Arc mode (expected MU, degrees)
TBI/TSET mode
PDD or TMR and OAF constancy
TBI/TSET output calibration
TBI/TSET accessories



- Doesn't have to be measured in water, but you already have your tank setup
- For a 10x10 field size, take measurements for monitor units of 500, 250, 100, 50, 25, 10, and 3 for each photon and electron beam energy
- Check to see if the dose is linear




TG-142: Annual Tests

X-ray/electron output calibration (TG-51)
Spot check of field size dependent output factors for x ray (two or more FSS)
Output factors for electron applicators (spot check of one applicator/energy)
X-ray beam quality (PDD₁₀ or TMR₁₀)
Electron beam quality (R_{50})
Physical wedge transmission factor constancy
X-ray monitor unit linearity (output constancy)
Electron monitor unit linearity (output constancy)
X-ray output constancy vs dose rate
X-ray output constancy vs gantry angle
Electron output constancy vs gantry angle
Electron and x-ray off-axis factor constancy vs gantry angle
Arc mode (expected MU, degrees)
TBI/TSET mode
PDD or TMR and OAF constancy
TBI/TSET output calibration
TBI/TSET accessories




- Doesn't have to be measured in water, but you already have your tank setup
- For a 10x10 field size, take measurements using 100 monitor units delivered at multiple dose rates for each photon and electron beam energy
- The readings should be the same




TG-142: Annual Tests

X-ray/electron output calibration (TG-51)
Spot check of field size dependent output factors for x ray (two or more FSS)
Output factors for electron applicators (spot check of one applicator/energy)
X-ray beam quality (PDD₁₀ or TMR₁₀)
Electron beam quality (R_{50})
Physical wedge transmission factor constancy
X-ray monitor unit linearity (output constancy)
Electron monitor unit linearity (output constancy)
X-ray output constancy vs dose rate
X-ray output constancy vs gantry angle
Electron output constancy vs gantry angle
Electron and x-ray off-axis factor constancy vs gantry angle
Arc mode (expected MU, degrees)
TBI/TSET mode
PDD or TMR and OAF constancy
TBI/TSET output calibration
TBI/TSET accessories



TG-142: Annual Tests


X-ray/electron output calibration (TG-51)
 Spot check of field size dependent output factors for x ray (two or more FSSs)
 Output factors for electron applicators (spot check of one applicator/energy)
 X-ray beam quality (PDD₁₀ or TMR₁₀)
Electron beam quality (R_{50})
 Physical wedge transmission factor constancy
 X-ray monitor unit linearity (output constancy)
 Electron monitor unit linearity (output constancy)
 X-ray output constancy vs dose rate
 X-ray output constancy vs gantry angle
 Electron output constancy vs gantry angle
 Electron and x-ray off-axis factor constancy vs gantry angle
 Arc mode (expected MU, degrees)
 TBI/TSET mode
 PDD or TMR and OAF constancy
 TBI/TSET output calibration
 TBI/TSET accessories



THE UNIVERSITY OF TENNESSEE

TG-142: Annual Tests


X-ray/electron output calibration (TG-51)
 Spot check of field size dependent output factors for x ray (two or more FSSs)
 Output factors for electron applicators (spot check of one applicator/energy)
 X-ray beam quality (PDD₁₀ or TMR₁₀)
Electron beam quality (R_{50})
 Physical wedge transmission factor constancy
 X-ray monitor unit linearity (output constancy)
 Electron monitor unit linearity (output constancy)
 X-ray output constancy vs dose rate
 X-ray output constancy vs gantry angle
 Electron output constancy vs gantry angle
 Electron and x-ray off-axis factor constancy vs gantry angle
 Arc mode (expected MU, degrees)
 TBI/TSET mode
PDD or TMR and OAF constancy
 TBI/TSET output calibration
 TBI/TSET accessories



THE UNIVERSITY OF TENNESSEE


TG-142: Annual Tests

X-ray/electron output calibration (TG-51)
 Spot check of field size dependent output factors for x ray (two or more FSS)
 Output factors for electron applicators (spot check of one applicator/energy)
 X-ray beam quality (PDD_{10} or TMR_{10}^{20})
 Electron beam quality (R_{50})
 Physical wedge transmission factor constancy
 X-ray monitor unit linearity (output constancy)
 Electron monitor unit linearity (output constancy)
 X-ray output constancy vs dose rate
 X-ray output constancy vs gantry angle
 Electron output constancy vs gantry angle
 Electron and x-ray off-axis factor constancy vs gantry angle
 Arc mode (expected MU, degrees)
 TBI/TSET mode
PDD or TMR and OAF constancy
 TBI/TSET output calibration
 TBI/TSET accessories




TG-142: Annual Tests

X-ray/electron output calibration (TG-51)
 Spot check of field size dependent output factors for x ray (two or more FSS)
 Output factors for electron applicators (spot check of one applicator/energy)
 X-ray beam quality (PDD_{10} or TMR_{10}^{20})
 Electron beam quality (R_{50})
 Physical wedge transmission factor constancy
 X-ray monitor unit linearity (output constancy)
 Electron monitor unit linearity (output constancy)
 X-ray output constancy vs dose rate
 X-ray output constancy vs gantry angle
Electron output constancy vs. gantry angle
 Electron and x-ray off-axis factor constancy vs gantry angle
 Arc mode (expected MU, degrees)
 TBI/TSET mode
PDD or TMR and OAF constancy
 TBI/TSET output calibration
 TBI/TSET accessories



10x10 with buildup cap



- For a 10x10 field size, take measurements using 100 monitor units delivered at gantry angles of 0, 90, 180, and 270 degrees for each photon beam energy
- Readings should be consistent with 10x10 baseline



TG-142: Annual Tests

X-ray/electron output calibration (TG-51)

Spot check of field size dependent
output factors for x ray
(two or more FSS)

Output factors for electron applicators
(spot check of one applicator/energy)

X-ray beam quality (PDD_{10} or TMR_{10}^{20})

Electron beam quality (R_{50})

Physical wedge transmission
factor constancy

X-ray monitor unit linearity
(output constancy)

Electron monitor unit linearity
(output constancy)

X-ray output constancy vs dose rate

X-ray output constancy vs gantry angle

Electron output constancy vs
gantry angle

Electron and x-ray off-axis factor
constancy vs gantry angle

Arc mode (expected MU, degrees)

TBI/TSET mode

PDD or TMR and OAF constancy

TBI/TSET output calibration

TBI/TSET accessories



- For a 10x10 field size, take measurements using 100 monitor units delivered at gantry angles of 0, 90, 180, and 270 degrees for each photon beam energy
- Readings should be consistent with 10x10 baseline



TG-142: Annual Tests

X-ray/electron output calibration (TG-51)

Spot check of field size dependent
output factors for x ray
(two or more FSS)

Output factors for electron applicators
(spot check of one applicator/energy)

X-ray beam quality (PDD_{10} or TMR_{10}^{20})

Electron beam quality (R_{50})

Physical wedge transmission
factor constancy

X-ray monitor unit linearity
(output constancy)

Electron monitor unit linearity
(output constancy)

X-ray output constancy vs dose rate

X-ray output constancy vs gantry angle

Electron output constancy vs
gantry angle

Electron and x-ray off-axis factor
constancy vs gantry angle


Arc mode (expected MU, degrees)

TBI/TSET mode


PDD or TMR and OAF constancy

TBI/TSET output calibration

TBI/TSET accessories



10x10 with
buildup cap



- For a 10x10 field size, take measurements using 100 monitor units delivered as 360 degree arcs and 180 degree arcs for each photon beam energy
- Readings should be consistent with 10x10 baseline



TG-142: Annual Tests

TABLE III. Annual.


Procedure	Machine-type tolerance		
	Non-IMRT	IMRT	SRS/SBRT
Mechanical			
Collimator rotation isocenter		± 1 mm from baseline	
Gantry rotation isocenter		± 1 mm from baseline	
Couch rotation isocenter		± 1 mm from baseline	
Electron applicator interlocks		Functional	
Coincidence of radiation and mechanical isocenter	± 2 mm from baseline	± 2 mm from baseline	± 1 mm from baseline
Table top sag		2 mm from baseline	
Table angle		1°	
Table travel maximum range movement in all directions		± 2 mm	
Stereotactic accessories, lockouts, etc.	NA	NA	Functional
Safety			
Follow manufacturer's test procedures		Functional	
Respiratory gating			
Beam energy constancy		2%	
Temporal accuracy of phase/amplitude gate on		100 ms of expected	
Calibration of surrogate for respiratory phase/amplitude		100 ms of expected	
Interlock testing		Functional	



THE UNIVERSITY OF TENNESSEE | OR

TG-142: Annual Tests

Mechanical	
Collimator rotation isocenter	
Gantry rotation isocenter	
Couch rotation isocenter	
Electron applicator interlocks	
Coincidence of radiation and mechanical isocenter	
Table top sag	
Table angle	
Table travel maximum range movement in all directions	
Stereotactic accessories, lockouts, etc.	
Safety	
Follow manufacturer's test procedures	
Respiratory gating	
Beam energy constancy	
Temporal accuracy of phase/amplitude gate on	
Calibration of surrogate for respiratory phase/amplitude	
Interlock testing	




- Using EPID or Film, overlapping fields with lower jaws at collimator angles of 0 and 180-degrees
- Use FWHM to examine the centering of the overlap



THE UNIVERSITY OF TENNESSEE | OR

TG-142: Annual Tests

- Mechanical**
- Collimator rotation isocenter
 - Gantry rotation isocenter
 - Couch rotation isocenter
 - Electron applicator interlocks
 - Coincidence of radiation and mechanical isocenter
 - Table top sag
 - Table angle
 - Table travel maximum range movement in all directions
 - Stereotactic accessories, lockouts, etc.
- Safety**
- Follow manufacturer's test procedures
- Respiratory gating**
- Beam energy constancy
 - Temporal accuracy of phase/amplitude gate on
 - Calibration of surrogate for respiratory phase/amplitude
 - Interlock testing




- Using EPID or Film, overlapping fields with upper jaws at collimator angles of 0 and 180-degrees
- Use FWHM to examine the centering of the overlap




THE UNIVERSITY OF TENNESSEE

TG-142: Annual Tests

- Mechanical**
- Collimator rotation isocenter
 - Gantry rotation isocenter
 - Couch rotation isocenter
 - Electron applicator interlocks
 - Coincidence of radiation and mechanical isocenter
 - Table top sag
 - Table angle
 - Table travel maximum range movement in all directions
 - Stereotactic accessories, lockouts, etc.
- Safety**
- Follow manufacturer's test procedures
- Respiratory gating**
- Beam energy constancy
 - Temporal accuracy of phase/amplitude gate on
 - Calibration of surrogate for respiratory phase/amplitude
 - Interlock testing



<input type="checkbox"/> Draw Beam Center on One View	<input type="checkbox"/> Draw Beam Center on Detail
<input type="checkbox"/> Beam Centers	<input type="checkbox"/> Beam Edges
<input type="checkbox"/> Beam Center Intersections	<input type="checkbox"/> Beam Intersections
<input type="checkbox"/> Farthest Beam Center Intersect Distance	<input type="checkbox"/> AND (AND) Edge Intersect Grid
<input type="checkbox"/> Minimum Target Circle (STAC)	<input type="checkbox"/> Maximum Hyperbola Distance Grid
<input type="checkbox"/> Max. Proj. Dist. Circle from User Defend Center	<input type="checkbox"/> Depth Planner
<input type="checkbox"/> Depth Per Line	<input type="checkbox"/> Beam Threshold Level -10



P1TG142TuxBrainFilm02072013Col StarComposite.nat	
<input type="checkbox"/> Beam Centers	<input type="checkbox"/> Beam Edges
<input type="checkbox"/> Beam Center Intersections	<input type="checkbox"/> Beam Intersections
<input type="checkbox"/> Beam Center Intersections	<input type="checkbox"/> Beam Center Intersections
<input type="checkbox"/> Farthest Distance Between Beam Center Intersections <0.5 mm	
Max. Dist. across Beam Edge Intersections = 3.5 mm	Mean Beam Width = 8.3mm
Max. Dist. across Beam Center Intersections = 3.5 mm	Mean Beam Width = 8.3mm
Minimum Proj. Distance to Beam from User Defined Center = 0.0mm	
Maximum Proj. Distance to Beam from User Defined Center = 0.16mm	

- Using EPID or Film, shoot 1-cm wide fields at collimator angles of 0, 45, 90, and 135 degrees
- Use FWHM to find center of “spokes” and look at variation



THE UNIVERSITY OF TENNESSEE

TG-142: Annual Tests

Mechanical

- Collimator rotation isocenter
- Gantry rotation isocenter
- Couch rotation isocenter
- Electron applicator interlocks
- Coincidence of radiation and mechanical isocenter
- Table top sag
- Table angle
- Table travel maximum range movement in all directions
- Stereotactic accessories, lockouts, etc.
- Safety**
- Follow manufacturer's test procedures
- Respiratory gating**
- Beam energy constancy
- Temporal accuracy of phase/amplitude gate on
- Calibration of surrogate for respiratory phase/amplitude
- Interlock testing

- Follow manufacturer's test procedures
- Respiratory gating
- Beam energy constancy
- Temporal accuracy of phase/amplitude gate on
- Calibration of surrogate for respiratory phase/amplitude
- Interlock testing

TG-142: Annual Tests

Mechanical

Collimator rotation isocenter
Gantry rotation isocenter
Couch rotation isocenter
Electron applicator interlocks
Coincidence of radiation and mechanical isocenter


Table top sag
Table angle
Table travel maximum range movement in all directions
Stereotactic accessories, lockouts, etc.

Safety

Follow manufacturer's test procedures

Respiratory gating

Beam energy constancy
Temporal accuracy of phase/amplitude gate on
Calibration of surrogate for respiratory phase/amplitude
Interlock testing



TG-142: Annual Tests

Mechanical

Collimator rotation isocenter
Gantry rotation isocenter
Couch rotation isocenter
Electron applicator interlocks
Coincidence of radiation and mechanical isocenter


Table top sag
Table angle
Table travel maximum range movement in all directions
Stereotactic accessories, lockouts, etc.

Safety

Follow manufacturer's test procedures

Respiratory gating

Beam energy constancy
Temporal accuracy of phase/amplitude gate on
Calibration of surrogate for respiratory phase/amplitude
Interlock testing



- On Elekta linacs, run the flexmap test



TG-142: Annual Tests

Mechanical

- Collimator rotation isocenter
- Gantry rotation isocenter
- Couch rotation isocenter
- Electron applicator interlocks
- Coincidence of radiation and mechanical isocenter**


- Table top sag
- Table angle
- Table travel maximum range movement in all directions
- Stereotactic accessories, lockouts, etc.

Safety

- Follow manufacturer's test procedures

Respiratory gating

- Beam energy constancy
- Temporal accuracy of phase/amplitude gate on
- Calibration of surrogate for respiratory phase/amplitude
- Interlock testing



- On Varian TrueBeam linacs, run the built-in MPC test



TG-142: Annual Tests

Mechanical

- Collimator rotation isocenter
- Gantry rotation isocenter
- Couch rotation isocenter
- Electron applicator interlocks**
- Coincidence of radiation and mechanical isocenter**


- Table top sag
- Table angle
- Table travel maximum range movement in all directions
- Stereotactic accessories, lockouts, etc.

Safety

- Follow manufacturer's test procedures

Respiratory gating

- Beam energy constancy
- Temporal accuracy of phase/amplitude gate on
- Calibration of surrogate for respiratory phase/amplitude
- Interlock testing



- Check to see that linac interlocks will clear only when the selected electron cone is inserted



TG-142: Annual Tests

Mechanical

Collimator rotation isocenter
Gantry rotation isocenter
Couch rotation isocenter
Electron applicator interlocks
Coincidence of radiation and mechanical isocenter
Table top sag
Table angle
Table travel maximum range movement in all directions
Stereotactic accessories, lockouts, etc.

Safety

Follow manufacturer's test procedures
Respiratory gating
Beam energy constancy
Temporal accuracy of phase/amplitude gate on
Calibration of surrogate for respiratory phase/amplitude
Interlock testing



- Check to see that linac interlocks will clear only when the selected stereotactic cone is inserted



TG-142: Annual Tests

Mechanical

Collimator rotation isocenter
Gantry rotation isocenter
Couch rotation isocenter
Electron applicator interlocks
Coincidence of radiation and mechanical isocenter
Table top sag
Table angle
Table travel maximum range movement in all directions
Stereotactic accessories, lockouts, etc.

Safety

Follow manufacturer's test procedures
Respiratory gating
Beam energy constancy
Temporal accuracy of phase/amplitude gate on
Calibration of surrogate for respiratory phase/amplitude
Interlock testing



- Place weight on table and check for sag
- Move table in all directions
- Rotate table and observe for walkout



TG-142: Annual Tests

Mechanical



Collimator rotation isocenter
 Gantry rotation isocenter
 Couch rotation isocenter
 Electron applicator interlocks
 Coincidence of radiation and mechanical isocenter
 Table top sag
 Table angle
 Table travel maximum range movement in all directions
 Stereotactic accessories, lockouts, etc.

Safety

Follow manufacturer's test procedures

Respiratory gating

Beam energy constancy
 Temporal accuracy of phase/amplitude gate on
 Calibration of surrogate for respiratory phase/amplitude
 Interlock testing



- Using a motion phantom of your choice to trigger the beam on and off
- Observe that the beam is turning on and off at the correct time
- Take measurements with an ionization chamber setup SAD at depths of 10 and 20-cm with the beam gating on and off for each energy
- Compare these readings to non-gated readings



TG-142: Dynamic/Virtual Wedges

TABLE IV. Dynamic/universal/virtual wedges.

Frequency	Procedure	Tolerance		
		Dynamic	Universal	Virtual
Daily	Morning check-out run for one angle		Functional	
Monthly	Wedge factor for all energies	C.A. axis 45° or 60° WF (within 2%) ^a	C.A. axis 45° or 60° WF (within 2%) ^a	5% from unity, otherwise 2%
Annual	Check of wedge angle for 60°, full field and spot check for intermediate angle, field size		Check of off-center ratios @ 80% field width @ 10 cm to be within 2%	

^aRecommendation to check 45° if angles other than 60° are used.




TG-142: Dynamic/Virtual Wedges


Daily Morning check-out run for one angle
Monthly Wedge factor for all energies

Annual Check of wedge angle for 60°, full field
and spot check for intermediate angle, field size


Physical



Virtual/Internal



Dynamic




TG-142: Dynamic/Virtual Wedges

Daily Morning check-out run for one angle
Monthly Wedge factor for all energies

Annual Check of wedge angle for 60°, full field
and spot check for intermediate angle, field size



- Include a 60-degree wedge angle on the morning QA routine
- Compare to baseline




TG-142: Dynamic/Virtual Wedges

Daily Morning check-out run for one angle
Monthly Wedge factor for all energies

Annual Check of wedge angle for 60°, full field
and spot check for intermediate angle, field size

- Include a 60-degree wedge angle in the monthly TG-51 test for each beam energy
- Use the open 10x10 to calculate the wedge factor and compare against baseline



THOMPSON CANCER SURVIVAL CENTER


THE UNIVERSITY OF TENNESSEE

TG-142: Dynamic/Virtual Wedges

Daily Morning check-out run for one angle
Monthly Wedge factor for all energies

Annual Check of wedge angle for 60°, full field
and spot check for intermediate angle, field size

- Using a diode array, chamber array, or EPID, measure a 60 degree and a 30 degree wedge angle profile for a 20x20 field size
- Compare against baseline



THOMPSON CANCER SURVIVAL CENTER

THE UNIVERSITY OF TENNESSEE

TG-142: MLC

TABLE V. Multileaf collimation (with differentiation of IMRT vs non-IMRT machines).

Procedure	Weekly (IMRT machines)	Tolerance
Qualitative test (i.e., matched segments, aka "picket fence")		Visual inspection for discernable deviations such as an increase in interleaf transmission
Setting vs radiation field for two patterns (non-IMRT)	Monthly	2 mm
Backup diaphragm settings (Elekta only)		2 mm
Travel speed (IMRT)		Loss of leaf speed >0.5 cm/s
Leaf position accuracy (IMRT)		1 mm for leaf positions of an IMRT field for four cardinal gantry angles. (<i>Picket fence</i> test may be used, test depends on clinical planning-segment size)
MLC transmission (average of leaf and interleaf transmission), all energies	Annually	±0.5% from baseline
Leaf position repeatability		±1.0 mm
MLC spoke shot		≤1.0 mm radius
Coincidence of light field and x-ray field (all energies)		±2.0 mm
Segmental IMRT (step and shoot) test		<0.35 cm max. error RMS, 95% of error counts <0.35 cm
Moving window IMRT (four cardinal gantry angles)		<0.35 cm max. error RMS, 95% of error counts <0.35 cm

We will cover these tests in a separate lecture



TG-142: Daily Imaging

TABLE VI. Imaging.

Procedure	Application-type tolerance	
	non-SRS/SBRT	SRS/SBRT
Daily ^a		
Planar kV and MV (EPID) imaging		
Collision interlocks	Functional	Functional
Positioning/repositioning	≤2 mm	≤1 mm
Imaging and treatment coordinate coincidence (single gantry angle)	≤2 mm	≤1 mm
Cone-beam CT (kV and MV)		
Collision interlocks	Functional	Functional
Imaging and treatment coordinate coincidence	≤2 mm	≤1 mm
Positioning/repositioning	≤1 mm	≤1 mm



TG-142: Daily Imaging

Planar kV and MV (EPID) imaging

Collision interlocks

Positioning/repositioning

Imaging and treatment coordinate coincidence
(single gantry angle)

Cone-beam CT (kV and MV)

Collision interlocks

Imaging and treatment coordinate coincidence

Positioning/repositioning



THOMPSON CANCER
SURVIVAL CENTER

THE UNIVERSITY OF TENNESSEE | OR

TG-142: Daily Imaging

Planar kV and MV (EPID) imaging

Collision interlocks

Positioning/repositioning



Imaging and treatment coordinate coincidence
(single gantry angle)

Cone-beam CT (kV and MV)

Collision interlocks

Imaging and treatment coordinate coincidence

Positioning/repositioning



THOMPSON CANCER
SURVIVAL CENTER

THE UNIVERSITY OF TENNESSEE | OR

TG-142: Monthly Imaging

TABLE VI. Imaging.

Procedure	Application-type tolerance	
	non-SRS/SBRT	SRS/SBRT
Monthly		
Planar MV imaging (EPID)		
Imaging and treatment coordinate coincidence (four cardinal angles)	≤2 mm	≤1 mm
Scaling ^b	≤2 mm	≤2 mm
Spatial resolution	Baseline ^c	Baseline
Contrast	Baseline	Baseline
Uniformity and noise	Baseline	Baseline
Planar kV imaging ^d		
Imaging and treatment coordinate coincidence (four cardinal angles)	≤2 mm	≤1 mm
Scaling	≤2 mm	≤1 mm
Spatial resolution	Baseline	Baseline
Contrast	Baseline	Baseline
Uniformity and noise	Baseline	Baseline
Cone-beam CT (kV and MV)		
Geometric distortion	≤2 mm	≤1 mm
Spatial resolution	Baseline	Baseline
Contrast	Baseline	Baseline
HU constancy	Baseline	Baseline
Uniformity and noise	Baseline	Baseline



TG-142: Monthly Imaging

Planar MV imaging (EPID)





- Imaging and treatment coordinate coincidence
(four cardinal angles)
- Scaling^b
- Spatial resolution
- Contrast
- Uniformity and noise

Planar kV imaging^d

- Imaging and treatment coordinate coincidence
(four cardinal angles)
- Scaling
- Spatial resolution
- Contrast
- Uniformity and noise

Cone-beam CT (kV and MV)

- Geometric distortion
- Spatial resolution
- Contrast
- HU constancy
- Uniformity and noise



TG-142: Monthly Imaging

Planar MV imaging (EPID)


Imaging and treatment coordinate coincidence
(four cardinal angles)
Scaling^b
Spatial resolution
Contrast
Uniformity and noise

Planar kV imaging^d

Imaging and treatment coordinate coincidence
(four cardinal angles)
Scaling
Spatial resolution
Contrast
Uniformity and noise

Cone-beam CT (kV and MV)

Geometric distortion
Spatial resolution
Contrast
HU constancy
Uniformity and noise



Leeds TOR 18FG Phantom



THE UNIVERSITY OF TENNESSEE | OR

TG-142: Monthly Imaging

Planar MV imaging (EPID)

Imaging and treatment coordinate coincidence
(four cardinal angles)
Scaling^b
Spatial resolution
Contrast
Uniformity and noise


Planar kV imaging^d

Imaging and treatment coordinate coincidence
(four cardinal angles)
Scaling
Spatial resolution
Contrast
Uniformity and noise

Cone-beam CT (kV and MV)

Geometric distortion
Spatial resolution
Contrast
HU constancy
Uniformity and noise

Low Contrast Resolution



THE UNIVERSITY OF TENNESSEE | OR

TG-142: Monthly Imaging

Planar MV imaging (EPID)


Imaging and treatment coordinate coincidence
(four cardinal angles)
Scaling^b
Spatial resolution
Contrast
Uniformity and noise

Planar kV imaging^d

Imaging and treatment coordinate coincidence
(four cardinal angles)
Scaling
Spatial resolution
Contrast
Uniformity and noise

Cone-beam CT (kV and MV)

Geometric distortion
Spatial resolution
Contrast
HU constancy
Uniformity and noise



TG-142: Monthly Imaging

Planar MV imaging (EPID)

Imaging and treatment coordinate coincidence
(four cardinal angles)
Scaling^b
Spatial resolution
Contrast
Uniformity and noise


Planar kV imaging^d

Imaging and treatment coordinate coincidence
(four cardinal angles)
Scaling
Spatial resolution
Contrast
Uniformity and noise

Cone-beam CT (kV and MV)

Geometric distortion
Spatial resolution
Contrast
HU constancy
Uniformity and noise

High Contrast Resolution



TG-142: Monthly Imaging

Planar MV imaging (EPID)

Imaging and treatment coordinate coincidence
(four cardinal angles)
Scaling^b
Spatial resolution
Contrast
Uniformity and noise



Planar kV imaging^d

Imaging and treatment coordinate coincidence
(four cardinal angles)
Scaling
Spatial resolution
Contrast
Uniformity and noise

Cone-beam CT (kV and MV)

Geometric distortion
Spatial resolution
Contrast
HU constancy
Uniformity and noise

Low Contrast Resolution



TG-142: Monthly Imaging

Planar MV imaging (EPID)

Imaging and treatment coordinate coincidence
(four cardinal angles)
Scaling^b
Spatial resolution
Contrast
Uniformity and noise



Planar kV imaging^d

Imaging and treatment coordinate coincidence
(four cardinal angles)
Scaling
Spatial resolution
Contrast
Uniformity and noise

Cone-beam CT (kV and MV)

Geometric distortion
Spatial resolution
Contrast
HU constancy
Uniformity and noise

Uniformity and Noise



$$\text{Spatial Uniformity} = (\text{CT}\#\text{max} - \text{CT}\#\text{min})/10$$
$$\text{Noise} = (\sigma * \text{CS} * 100) / \mu\text{w}$$



Questions



THOMPSON CANCER
SURVIVAL CENTER

THE UNIVERSITY OF TENNESSEE | UT

TG-142 vs TG-40

Monthly		
	TG-40	TG-142 changes
Dosimetry		
x-ray central axis dosimetry parameter (PDD, TAR) constancy	2%	Removed
Electron central axis dosimetry parameter constancy (PDD)	2mm	2%/2mm
x-ray beam flatness constancy	2%	<i>Replaced with 1% constancy of profile</i>
Electron beam flatness constancy	3%	
x-ray and electron symmetry	3%	
Interlock Checks		
Emergency Off	Functional	Removed
Wedge, "cone"	Functional	
Mechanical		
Light/radiation field coincidence	2 mm or 1%/side	<i>Only if clinical setups performed</i>
Field size indicators	2mm	1mm/side
Cross-hair centering	2mm	1mm
Treatment couch position indicators	2 mm/l deg	Tighter for SRS/SBRT



THOMPSON CANCER
SURVIVAL CENTER

THE UNIVERSITY OF TENNESSEE | UT

TG-142 vs TG-40

Annual		
	TG-40	TG-142 changes
Dosimetry		
x-ray/electron output calibration constancy	2%	+,- 1%
Field size dependence of x-ray output constancy	2%	Spot check 2 or more FS
Output factor constancy for electron applicators	2%	Combination of energies & applicators (i.e., Varian -5)
Central axis parameter constancy (PDD, TAR)	2%	+ 1% for TG-51 purpose
Off-axis factor constancy	2%	+ 1% from baseline
Wedge transmission factor constancy	2%	Unique for wedge systems
Monitor chamber linearity	2%	Unique for delivery system
x-ray output constancy vs gantry angle	1%	
Electron output constancy vs gantry angle	2%	+ 1% from baseline
Off-axis factor constancy vs gantry angle	2%	



THE UNIVERSITY OF TENNESSEE | OR

TG-142: Annual Imaging

TABLE VI. Imaging.

Procedure	Application-type tolerance	
	non-SRS/SBRT	SRS/SBRT
Annual (A)		
Planar MV imaging (EPID)		
Full range of travel SDD	±5 mm	±5 mm
Imaging dose ^c	Baseline	Baseline
Planar KV imaging		
Beam quality/energy	Baseline	Baseline
Imaging dose	Baseline	Baseline
Cone-beam CT (kV and MV)		
Imaging dose	Baseline	Baseline



THE UNIVERSITY OF TENNESSEE | OR

Questions



THOMPSON CANCER
SURVIVAL CENTER

THE UNIVERSITY OF TENNESSEE | The logo for The University of Tennessee, featuring the letters "ut" in orange and white.